

Brescia, 18 January 2019

DECLARATION ON eSWAB TIP COLOR


Copan Italia S.p.A. declares that:

The flocked swabs included in the eSwab product are coated to assist in the preservation of bacteria. This coating, together with a technical improvement in the air flow of the oven used to dry the tips during the production of the swabs, can cause the fibers of the flocked swab to naturally turn from light yellow to brownish color. However, the tip color does not in any way affect the performance of the swab.

Yours Faithfully,



Alessia Tomasi
Copan Group Chief Quality Officer



copanitalia.com